

RAPID AND AFFORDABLE POINT-OF-CARE MOLECULAR DIAGNOSTICS



# **Lodestar DX**

## Molecular UTI testing at the point-of-care

- 98% sensitivity and 91% specificity for *E.coli*
- 88% overall sensitivity and 87% specificity
- Results in 35 minutes
- No extraction step required
- Test for six pathogens simultaneously
- Small and portable fits into the busiest lab, clinic or office
- Test the patient sample at any time on the day of collection





### Lodestar DX

- Isothermal LAMP technology provides molecular accuracy comparable to microbiological culture analysis
- Small footprint: 140 mm (w) x 200 mm (d) x 90 mm (h)
- Lightweight: 720 g
- Simple sampling steps with minimal preparation
- Green light indicates a negative result; red light indicates a positive result
- Operates on 240-12V power
- Connect to LIMS via Bluetooth or by USB
- No analyser maintenance or calibration is required
- 12 month warranty from date of installation

#### **UTI test panel**

- Test for the six most common pathogens that cause UTI from one sample
- Detects multiple pathogens in a mixed growth sample
- Test strip includes positive and negative controls
- Only 10 µL of urine sample is required
- Allows informed prescribing decisions to be made

### **Clinical performance evaluation**

Assay	Sensitivity %	Specificity %
E. coli	98	91
Enterococcus/S. aureus	71	95
Staph. saprophyticus	100	98
Proteus mirabilis	86	96
Pseudomonas aeruginos	a 100	96
Klebsiella pneumoniae	100	98

# Test urine for six pathogens quickly and easily



Add urine to the diluent tube



Mix into all six assay tubes and positive control



Place test panel in Lodestar DXResults in 35 minutes

Llusern Scientific Cardiff Edge Business Park Longwood Drive, Cardiff CF14 7YU United Kingdom sales@llusern.co.uk

www.llusern.co.uk

HUTI-EN-PS-03